

Phonak Audeo Fit™

Product information



Our latest innovation in the Paradise line, Audeo Fit has a heart rate sensor that can track the body's health data which can be viewed with the myPhonak app. By tracking heart rate, steps and activity levels, these hearing aids can promote healthy habits and create a greater sense of well-being. Suitable for mild to severe hearing losses and powered by the innovative technology of the Phonak Paradise platform, Phonak Audeo Fit hearing aids are the solution for those clients that want to engage with their health data.

For more information visit www.phonakpro.com

A Sonova brand

PHONAK
life is on

Product description



Audéo P-R Fit
Audéo P-R Fit Trial
P90/P70

Multi-function button		•
Bluetooth®		•
Hydrophobic coating		•
IP rating		IP68 ¹
Dimensions	L x W x D	30.6 x 12.3 x 8.1 mm (1.21 x 0.49 x 0.32")
	Weight	2.23 g (0.078 oz)
Charging times ²	0% to 100%	3 h
	0% to 80%	1.5 h
	30 min. fitting	10 min.
M Receiver		
Maximum Power Output (dB SPL)		114
Max. gain (dB)	2 cc coupler	51
Frequency range (Hz)	2 cc coupler	<100 - >8000
P Receiver		
Maximum Power Output (dB SPL)		122
Max. gain (dB)	2 cc coupler	59
Frequency range (Hz)	2 cc coupler	<100 - 6300

Housing colors



P1
Sand Beige



P4
Chestnut



P6
Silver Gray

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 meter of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard.
² Actual charging time depends on the remaining battery life, but it will not exceed 3 hours.

Performance levels

	Premium	Advanced	Features	Premium	Advanced
AutoSense OS™ 4.0	Premium	Advanced	UltraZoom	•	•
Comfort in Echo	•		Dynamic Noise Cancellation	•	
Speech in Car	•		Speech Enhancer	•	
Speech in Loud Noise	•		EchoBlock	•	
Music	•	•	Tap control		
Comfort in Noise	•	•	- Siri/Google Assistant	•	•
Speech in Noise	•	•	- Pause/resume streaming	•	•
Calm Situation	•	•	- Accept/end phone call	•	•
Media speech + mic	•	•	StereoZoom	•	•
Media music + mic	•	•	WindBlock	•	•
Additional programs			SoundRelax	•	•
Max. additional programs	4	4	Motion Sensor Hearing	•	•
Comfort in Echo	•		DuoPhone	•	•
Speech in Loud Noise	•	•	Real Ear Sound	•	•
Speech in 360°	•	•	Roger and Directional setting	•	•
Speech in Noise	•	•	SoundRecover2	•	•
Calm Situation	•	•	WhistleBlock	•	•
Comfort in Noise	•	•	NoiseBlock	•	•
Music	•	•	RogerDirect™	•	•
Acoustic phone	•	•	User Preference Tuning	•	•
Custom program	•	•	QuickSync	•	•
Streaming programs			AOV	•	•
RogerDirect + mic	•	•	Tinnitus Balance	•	•
PartnerMic + mic	•	•	auto Acclimatization	•	•
Phone call + mic	•	•	Environmental balance	•	•
			Fine tuning channels	20	20

Health features

	Audéo Fit	Paradise rechargeable RICs
Separate wearing time for each hearing aid	•	•
Average combined wearing time for bilateral hearing aids	•	•
Average wearing time in different sound environments	•	•
Steps	•	•
Activity levels	•	•
Distance walked or ran	•	
Heart rate	•	
Optional goal setting for wearing time	•	•
Optional goal setting for steps	•	•

The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc.
IOS is a trademark of Cisco, Inc.
Android is a trademark of Google, Inc.

Fitting

Audéo P-R Fit Audéo P-R Fit Trial

Phonak Target	7.3 or higher
Noahlink Wireless	•

Acoustic coupling*

Length

Available in four different lengths (0, 1, 2, 3), left and right

M Sensor Receiver

•

P Sensor Receiver

•

Dome

Cap dome

•

Open dome Available in three sizes (S, M, L)

•

Vented dome clear Available in three sizes (S, M, L)

•

cShell


Acryl

•


•

Fitting ranges

* Only the 8 pin Sensor Receiver and Sensor cShell works with Audeo Fit.



Mild to severe hearing loss.



Mild to severe hearing loss.

Phonak Audéo Fit Trial™

The new Phonak Audéo P Fit Trial hearing aids can be quickly set to the required performance level in order to suit your recommendation. Both levels of technology in one Phonak Audéo P Fit Trial hearing aid provides you with many options for providing the hearing solution at the first appointment or bridging the gap during repairs.

Phonak Audéo Fit charging options



Phonak Charger Case Combi

The Phonak Charger Case Combi is a charger, drying kit and protective hard case all in one. Additionally, it includes a cleaning tool in the robust and splash-proof housing.



Phonak Power Pack

For charging rechargeable hearing aids on the go, the Phonak Power Pack can be easily attached to the bottom of the Phonak Charger Case Combi.

Phonak Audéo Fit Bluetooth® versions

Direct calls/audio streaming

Bluetooth® 4.2 wireless technology and most older Bluetooth phones

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc.

Wireless communication portfolio

	Audéo P-R Fit Audéo P-R Fit Trial P90/P70
Wireless accessory	
TV Connector	•
Phonak PartnerMic	•
Phonak RemoteControl	•
Roger	
RogerDirect™	•
Apps	
myPhonak	•

For more information please visit www.phonakpro.com/audeoparadise

For information on Roger microphones please visit www.phonakpro.com/roger

For more information please visit www.phonakpro.com/smartapps

life is on


At Phonak, we believe that well-hearing equates to well-being and is essential to living life to the fullest. For more than 70 years, we have remained passionate about creating a world where 'life is on' for everyone. Our innovative hearing solutions are designed for people of all ages and all degrees of hearing loss, to connect socially, thrive mentally and emotionally.

www.phonak.com



027-0666-02/V1.00/2022-01/na © 2022 Sonova AG All rights reserved



 Sonova AG · Laubisrütistrasse 28
CH-8712 Stäfa · Switzerland

